

Background

Endophthalmitis is the most challenging complication seen in ophthalmic practice. Culture of intraocular specimens like aqueous humour, vitreous humour and vitrectomy fluid, is still considered the gold standard in the diagnosis of endophthalmitis. However, the diagnostic yield of culture method is only 25-60% of the clinically diagnosed typical cases (Long C 2014)(Nicoara SD 2014).

Recent study showed that metagenomic next-generation sequencing (mNGS) has a higher positive rate of identifying pathogens in endophthalmitis than in culture (Zhu, et al. 2022). The positive rates of mNGS and culture were 88.89% (32/36) and 27.78% (10/36), respectively (Zhu, et al. 2022). *Staphylococcus epidermidis*, *Streptococcus pneumoniae*, and *Klebsiella pneumoniae* were the most pathogenic bacteria in traumatic, postoperative, and endogenous endophthalmitis, respectively.

Intraocular fluids for microbial testing in endophthalmitis cases are usually limited in sample amounts, which render the sample processing rather challenging. Devin™ Fractionation Filter utilizes Zwitterionic interface Ultra-Self-assemble Coating Technology (ZISC Technology) to specifically binds to human nucleated cells in the human biospecimen while the 15-20 µm membrane pores allow other components including microorganisms to pass through. The objective of this study is to test if the sample of small volume can be processed with PaRTI-Seq workflow with a modified protocol.

Material and Methods

A sample of vitreous humour (labelled VH; from infected eye) and an eye swab (labelled UI; from uninfected eye) were collected. The eye swab was intended as an “uninfected” control, representative of the microflora of the patient’s eye surface, for comparison with the intraocular VH sample. The samples were dispatched immediately after collection and transported to Micronbrane’s laboratory at temperature below 10°C.

The fluidic samples were processed with a modified protocol using Devin™ Fractionation Filter for samples of small volumes as illustrated in **Fig. 1**. First, the Devin™ filter was pre-wet with 500 µL of PBS. To the 50 µL VH sample received, 450 µL of PBS was added to make up to a volume of 500µL. The eye swab was immersed in 500 µL PBS. Both samples were then filtered through Devin™ filter using 3 mL syringe. Another 500 µL PBS wash was applied to filter to collect filtrations of both steps.

The final volume of the sample after filtration was about 1.0ml. The samples were then subject to centrifugation at 16,000g for 15 mins. DNA was purified from the microbial pellet and constructed into NGS library using Unison™ ultralow DNA NGS library kit for Illumina platforms*.

A non-template control (NTC) was prepared in parallel from the time of DNA extraction. Approximately 20 million NGS reads were obtained for both samples and analysed by in-house bioinformatic pipeline. The diagnostic performance was compared between mNGS and culture results from the same sample.




Fig. 1 PaRTI-Seq™ workflow with modified protocol for the liquid biopsy samples of small volume (< 500 µL) described in detail

*All the reagents and consumables required for DNA extraction and library construction for the PaRTI-Seq™ workflow are supplied by Micronbrane.

Result

NGS sequences of 100bp were used for the analysis in this application. The pipeline carries out human DNA subtraction and reports the number of microbial reads (Micro_Reads) mapped to the ~1400 microorganisms in the database. The pipeline also reports the percentage of Microbial reads (Micro_%) of each species among all the microorganisms reads.

The results of primary analysis of the samples and NTC are shown in **Table 1** below. The Micro_%(SPL): Micro_%(NTC) ratio of each species is computed for both samples. The top 10 species with the highest ratio is listed in **Table 1**. A threshold of 10 was arbitrarily set as the threshold for a valid positive call. As a result, the top positive call was *Klebsiella pneumoniae*, which was concordant with the identification of the same organism by microbial culture.

A plot of the Micro_%(SPL): Micro_%(NTC) ratio in the VH sample (**Fig. 2**) shows that a true positive is likely to have ratios of up to 1000 and can be easily identified from background noise. Interestingly, other *Klebsiella* spp such as *K. quasipneumoniae* and *K. varicola* were also called alongside *K. pneumoniae*. The alignment to these non-*K. pneumoniae* *Klebsiella* spp are probably due to sequences that are conserved across the *Klebsiella* genus and sequences were aligned using only 100 bp of read length.




Fig 2. Plot of Micro%Sample:Micro%NTC ratio in VH and UI

However, as (Schlaberg 2017) pointed out, it is not possible to distinguish between partial coverage due to low bacterial genome copies and that due to conserved sequences. Another possibility is that these represent genuine species heterogeneity in a mixed infection that is easily missed by traditional culture.

Conclusion

The Endophthalmitis Vitrectomy Study (EVS) in 1995 showed that nearly 30% of vitreous culture specimens were culture negative (Group 1995). Subsequent studies have had more variable success, with culture-positive rates ranging from 38 to 64% from post trauma endophthalmitis for example (Long C 2014) (Nicoara SD 2014). The actual clinical performance of mNGS for eye infection is not known but is likely to have a higher diagnostic yield (ie. positive rates) as compared to traditional methods, like studies with other sample types have shown.

This study highlighted that small volume sample can be processed with PaRTI-Seq™ workflow with a modified protocol. In conclusion, we demonstrated Devin™ filter can be used to process biospecimen of small volume as low as 50 µL. The following DNA extraction and NGS library preparation reagents in PaRTI-Seq™ workflow were all compatible with low sample inputs such as vitreous humour samples. By extension, PaRTI-Seq™ is likely to be compatible with other scarce samples such as CSF.

Table 1: Primary analysis results of VH and UI samples compared with NTC. Top 10 species with the highest Micro_%(SPL):Micro_%(NTC) ratio (largest to smallest) is listed. A threshold of 10 was arbitrarily set as the threshold for a valid positive call. Positive calls are highlighted in yellow

Organism	Sample VH			Sample UI			NTC	
	Micro_Reads	Micro_%	Micro_%(SPL): Micro_%(NTC) ratio	Micro_Reads	Micro_%	Micro_%(SPL): Micro_%(NTC) ratio	Micro_Reads	Micro_%
<i>Klebsiella pneumoniae</i>	2689	3.9795	1047.24	13	0.0168	4.42105	116	0.0038
<i>Klebsiella quasipneumoniae</i>	354	0.5239	141.595	13	0.0168	4.42105	112	0.0037
<i>Klebsiella variicola</i>	399	0.5905	74.7468	19	0.0246	3.11392	239	0.0079
<i>Klebsiella aerogenes</i>	127	0.1879	21.5977	21	0.0272	3.12644	264	0.0087
<i>Citrobacter amalonaticus</i>	34	0.0503	6.36709	37	0.0479	6.06329	240	0.0079
<i>Malassezia globosa</i>	31	0.0459	5.53012	31	0.0401	4.83133	253	0.0083
<i>Burkholderia contaminans</i>	10917	16.1561	4.70091	11311	14.629	4.25658	104592	3.4368
<i>Mycobacteroides saopaulense</i>	46	0.0681	4.57047	58	0.075	5.03356	453	0.0149
<i>Mycobacteroides franklinii</i>	368	0.5446	4.41687	396	0.5122	4.1541	3752	0.1233
<i>Pseudomonas fulva</i>	30	0.0444	4.18868	51	0.066	6.22642	322	0.0106

References

- Afshinnekoo, et al. "Precision Metagenomics: Rapid Metagenomic Analyses for Infectious Disease Diagnostics and Public Health Surveillance." J Biomol Tech 28 (1), 40–45. doi: 10.7171/jbt.17-2801-007, 2017: 40-45.
 Group, Endophthalmitis Study. "Results of the Endophthalmitis Vitrectomy Study." Arch Ophthalmol 1995; 113:1479–1496., 1995: 1479-1496.
 Long C, et al. "Causative organisms of posttraumatic endophthalmitis: a 20-year retrospective study." BMC Ophthalmol 2014; 14:34., 2014: 34.
 Nicoara SD, et al. "Outcome and prognostic factors for traumatic endophthalmitis over a 5-year period." J Ophthalmol 2014; 2014: 747015, 2014: 747015.
 Simner, et al. "Understanding the Promises and Hurdles of Metagenomic Next-Generation Sequencing as a Diagnostic Tool for Infectious Diseases." Clin Infect Dis 66 (5), 778–788. doi: 10.1093/cid/cix8, 2018: 778-788.
 Zhu, Junfeng BS, et al. "Metagenomic next-generation sequencing detects pathogens in endophthalmitis patients." Retina: May 2022 - Volume 42 - Issue 5 - p 992-1000, 2022: 992-1000.

For Research Use Only

Micronbrane Medical | 22F., No. 99, Xinpu 6th St., Taoyuan Dist., Taoyuan City 33044, Taiwan | info@micronbrane.com

Tel: 03-316-6428 | Fax 03-316-6408